# SHIMADZU

### The rapid detection of botanical natural toxins using probe ESI unit combined with quadrupole time-of-flight mass spectrometer

<u>Tetsuo lida<sup>1</sup>, Kaoru Nakagawa<sup>1</sup>, Manami Kobayashi<sup>2</sup></u> <sup>1</sup> Shimadzu Corporation, Kyoto-city, Japan, <sup>2</sup> Shimadzu Corporation, Kawasaki-city, Japan

## 1. Introduction

Many animals and plants are known to have poisonous components (natural toxins) in their bodies. Food poisoning caused by animals and plants containing these natural toxins occurs every year, although the number of cases are not very high compared to bacterial food poisoning. However, this is extremely important from the viewpoint of food sanitation because some poisons have a high fatality rate. Therefore, rapid and effective screening methods are required to guarantee food safety. In this study, we report a method for rapid detection of natural plant toxins in Colchicum autumnale and Narcissus tazetta using DPiMS<sup>™</sup> QT equipped with probe electrospray ionization (PESI) technology. Additionally, we use a DPiMS QT mounted on the LCMS-9050, a newly designed quadrupole time-of-flight (Q-TOF) mass spectrometer that can rapidly achieve mass stability. This combination enables direct analysis and minimizes the time from sample preparation to accurate analysis.

## 2. Sample Preparation and Analysis Conditions

The leaves and bulbs of Narcissus tazetta and the bulbs of Colchicum autumnale were obtained from garden stores. Fig.1 shows the natural toxins contained in these plants. The sample was cut into approximately 5 mm squares with a thickness of 1 to 2 mm and set in a biological sample plate. 35  $\mu$ L of 50% ethanol-water solution (v/v) was added to the well of this plate. In the PESI-Q-TOF system (Shimadzu Corporation), a probe picks up sample from the sampling plate. At the same time, by applying a voltage to the probe, the sample adhering to the probe surface is ionized and directly introduced into the mass spectrometer (Fig. 2). Table 1 shows the analytical conditions for this analysis. Qualitative analysis was performed using the LabSolutions Insight Explore<sup>TM</sup> software (Shimadzu Corporation).



Fig. 1 Structural formula of alkaloids in Colchicum autumnale and Narcissus tazetta



### Analysis of natural toxins contained in Colchicum autumnale

Fig. 3 shows the mass spectrum of a *Colchicum autumnale* bulb section measured in the positive mode. By using a dedicated analytical software, component A (m/z 400.1756) was estimated to be  $C_{22}H_{25}NO_6$  with high mass accuracy of within 1 mDa compared to the theoretical mass value (Fig. 4). When an online search was conducted using the ChemSpider database for the estimated composition formula, colchicine was suggested as a top candidate compound. Fig. 5 shows the automatically detected match between the product ions predicted from the structural formula and the product ions observed in the MS/MS spectrum. These results strongly support that component A is colchicine. In the same manner, hexose-based saccharides, metabolites choline and arginine, and the natural toxin demecolcine were also detected.

Fig. 5 Automatic assignment of MS/MS Fragments

### Analysis of natural toxins contained in Narcissus tazetta

Fig. 6 shows the mass spectrum of a Narcissus tazetta bulb section measured in the positive mode. Accurate mass spectrometry was performed using analysis software, showing that saccharides composed of hexose, metabolites choline and arginine, and natural toxins tazettine, lycorine, and galanthamine, were detected. Figure 7 shows that MS/MS spectral pattern of lycorine obtained from the bulb measurement matches that of the standard. These components was also detected in a section of *Narcissus tazetta* leaf.

mass accuracy.

• The time required for pretreatment was approximately 5 minutes and the measurement time was 0.5 minutes, significantly reducing the time required for analysis compared to measurement by LC or LC/MS.

 The obtained data can be used for simple and rapid composition estimation and structural analysis of unknown components using the qualitative analysis software.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures. DPiMS and LabSolutions Insight Explore are rademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.