

The 4th International Electronic Conference on Foods , 15–30 Oct 2023



Evaluation of DNA extraction methods for PCR analysis of maize and sunflower oils

**Tamara Kutateladze¹, Kakha Bitskinashvili¹, Boris Vishnepolsky¹, Kakha Karchkhadze², Tata Ninidze^{1, 2},
Nelly Datukishvili^{1, 2}**

¹ **Ivane Beritashvili Center of Experimental Biomedicine, Tbilisi, Georgia**

² **School of Natural Sciences and Medicine, Iliava State University, Tbilisi, Georgia**

Challenges of the Modern Food Industry



The modern food industry has a big challenge dealing with the adulteration of different types of plant oils

Oil plants differ in the following indicators:

- Chemical and molecular composition
 - Toxic and allergenic properties
 - Nutritional value and safety
 - industrial value



- **International legislation requires food authentication and traceability at all stages of production**
- It is necessary to provide accurate information about the presence of food ingredients and allergens through labeling

A food monitoring system needs effective methods of detecting food ingredients!

Study Aim:

To develop efficient PCR technology for plant DNA traceability in oils

Objects:

Sunflower and Maize

Seeds

Flours

cold-pressed oils

refined oils

Sunflower samples



Sunflower and maize seeds and flours were purchased from a supermarket



Cold-pressed sunflower oil was provided by the company „AgroPro Ltd”

Maize samples






Refined and cold-pressed corn oils were purchased from a supermarket



Preparation of DNA samples of appropriate quality and sufficient quantity from oils remains a major drawback for successful PCR detection!


Our study examines several approaches to DNA enrichment and extraction



International Journal of Food Science & Technology
Short communication
Comparative evaluation of DNA extraction methods for food crops
Nelly Datskibidze,^{1,2*} Inga Gabriadze,¹ Tamara Kutateladze,¹ Marina Karseladze,¹ Boris Vishnepolsky¹
¹ Institute of Molecular Biology and Biotechnological Physics, 12, Gruzia St., Tbilisi 0100, Georgia
² Ilia Chavchavadze State University, 52 Chavchavadze Ave., Tbilisi 0179, Georgia



Genomic DNA extraction



Agarose gel electrophoresis of genomic DNAs


Oil DNA enrichment – Centrifugation of 30 ml cold-pressed sunflower oil at 18,000 g, 4 °C

DNA samples:


1. cold-pressed sunflower oil - extraction by modified CTAB method
2. cold-pressed sunflower oil - extraction by standard CTAB method
- 3 - 4. sunflower seeds - extraction by standard CTAB method
5. Maize flour - extraction by standard CTAB method

The genomic DNA band visible on the agarose gel was obtained from cold-pressed sunflower oil!

PCRs specific to 18S ribosomal RNA gene



Amplicon 18S-302



Amplicon 18S-140

Samples


- 1-2. cold-pressed sunflower oil, Olive Oil kit
- 3-4. cold-pressed sunflower oil, NucleoSpin Food kit
5. Sunflower seeds, plant mini kit
6. Water negative control

Results

- DNA enrichment by oil centrifugation at 18,000 g, 4 °C gave best results
- Olive oil kit and Nucleospin food kit gave amplifiable DNA from cold pressed sunflower oil

PCRs specific to 18S ribosomal RNA gene


Comparison standard and modified CTAB methods



Amplicon 18S-302

Samples:


1. Sunflower seeds, **standard CTAB method**
2. cold-pressed sunflower oil, **standard CTAB method**
3. cold-pressed sunflower oil, **modified CTAB method**
4. Water -negative control



Amplicons 18S-140 (1-5) and 18S – 167(6-10)

Samples:

- 1, 6. Sunflower seeds, **standard CTAB method**
- 2, 7. cold-pressed sunflower oil, **standard CTAB method**
- 3-4, 8-9. cold-pressed sunflower oil, **modified CTAB method**
- 5.10. Water negative control



Amplicons 18S-140 (1-4) and 18S – 167 (5-8)

Samples:


- 1, 5. maize seeds, standard CTAB method
- 2, 6. refined maize oil, modified CTAB method
- 3, 7. cold-pressed sunflower oil, modified CTAB method
- 4, 8. Water negative control

Results

- The modified CTAB method gave amplifiable DNA from oils, but standard CTAB method did not give genomic DNA from oils
- Standard CTAB method gave amplifiable DNA from Sunflower and maize seeds
- The expected amplicons were amplified in both seeds and oil samples by PCRs targeting 18S RNA gene
- 18S302 amplicon-specific PCR is a more efficient method than 18S167 and PCR 18S140 amplicon-specific PCRs

Sunflower detection in oils by PCR

Comparison of standard and modified CTAB methods




PCR amplicons
heli-77 (1-4); heli-104 (6-9); heli-160 (11-14)



PCR amplicons
heli-83 (1-4); heli-188 (6-9)



PCR amplicons
heli-162 (1-4); heli-188 (6-9)



PCR amplicons
Heli-160 (1-4), heli-162 (5-8);
heli-188 (9-12); heli-104 (13-16)



PCR amplicons
Heli-77 (1-4), heli-83 (5-8);
heli-160 (9-12)

Samples

- 1, 6, 11. Sunflower seeds, standard CTAB method
- 2, 7, 12. 2.cold-pressed sunflower oil, **standard CTAB method**
- 3, 8, 13. cold-pressed sunflower oil, **modified CTAB method**
4. Water negative control

Samples

- 1, 5, 9, 13. cold-pressed sunflower oil, **Olive Oil kit**
- 2, 6, 10, 14. cold-pressed sunflower oil, **NucleoSpin Food kit**
- 3, 7, 11, 15. Sunflower seeds, plant mini kit
- 4, 8, 12, 16. Water negative control

Results

- A modified CTAB method was found to be the best method for extracting amplifiable DNA from sunflower oil
- The standard CTAB method was found to be the best method for extracting amplifiable DNA from sunflower seeds
- DNA obtained by the NucleoSpin Food Mini Kit is useful for PCR detection of sunflower in oil
- DNA obtained with the Olive Oil DNA Isolation Kit failed to detect sunflower by PCR in the oil

Conclusions

- The oil type (cold-pressed, refined), DNA extraction and amplification methods are important for successful PCR analysis
- Centrifugation of the oils at 18,000 g, at 4 °C was the best method for DNA enrichment from the oils
- The modified CTAB method was found to be the best DNA extraction method for PCR analysis of sunflower and maize oils
- In addition, a PCR system specific for the 18S-302 amplicon of the 18S ribosomal RNA gene was identified as the best method for DNA traceability in oils



Acknowledgments

- This work was supported by Shota Rustaveli National Science Foundation of Georgia (SRNSFG) [Grant AR-22- 636]
- The co-funding organization of the project is "Biodiesel Georgia" LLC.

Thank You for Your Attention!



Ilia State University

Nelly Datukishvili

Kakha Karchkhadze

Tata Ninidze

Tamar Koberidze

*I. Beritashvili Center of
Experimental Biomedicine*

Kakha Bitskinashvili

Tamar Kutateladze

Boris Vishnepolsky

