


1ST INTERNATIONAL ELECTRONIC CONFERENCE ON CHEMICAL SENSORS AND ANALYTICAL CHEMISTRY


01-15 JULY 2021 | ONLINE




TROPOMYOSIN ANALYSIS IN FOODS USING AN ELECTROCHEMICAL IMMUNOSENSING APPROACH

Ricarda Torre, Maria Freitas, Estefanía Costa-Rama, Henri P.A. Nouws, Cristina Delerue-Matos





requimte | Laboratório Associado



LNC

Analytical features

Applicability to commercial food samples

Linear Range	2.5 – 20 ng mL ⁻¹
Correlation coefficient	0.990
Slope (m)	0.787 µA ng⁻¹ mL
LOD	1.7 ng mL ⁻¹
LOQ	5.7 ng mL ⁻¹
V _{x0}	8.8 %

Sample	Detected value
Shrimp	$80.42 \pm 2.7 \ \mu g \ g^{1}$
Shrimp sauce	$170.4 \pm 1.80~\text{ng g}^{-1}$
Crab paste	$21.6\pm4.13~\text{ng g}^{-1}$
Chicken paste	n.d.
n.d. not detected	

Final Remarks

- A simple immunosensor for tropomyosin analysis was developed;
- This immunoassay only takes 2h50 min, and it requires 40 µL of sample to perform the analysis;
- The sensor can determine tropomyosin in a concentration range between 2.5 and 20 ng mL⁻¹ with a good precision (coefficient of variation of 8.8 %);
- It was successfully applied to commercial food samples. Previous values of tropomyosin will be compared with traditional ELISA methods in further studies.

Funding

This research was funded by Portuguese national funds (FCT/MCTES, Fundação para a Ciência e a Tecnologia and Ministério da Ciência, Tecnologia e Ensino Superior) through the project UIDB/50006/2020. This work also received financial support from the European Union (FEDER funds through COMPETE POCI-01-0145-FEDER-030735) and National Funds (FCT) through the project PTDC/QUI-QAN/30735/2017 - TracAllerSens - Electrochemical sensors for the detection and quantification of trace amounts of allergens in food products.R.Torre thanks FCT and ESF (European Social Fund) through POCH (Programa Operacional Capital Humano) for her PhD grant (SFRH/BD/143753/2019).