

ULTRALIGHT, COMPACT, AND STRETCHABLE ELECTRONICS FOR CONTINUOUS AND HIGH-QUALITY CARDIAC ASSESSMENT IN HIGH USER ACTIVITIES

IEN Center for Human-Centric
Interfaces and Engineering (CHCIE)

Yun-Soung Kim, Nathan A. Rodeheaver, Musa Mahmood, Hojoong Kim, and Woon-Hong Yeo
George W. Woodruff School of Mechanical Engineering,
Wallace H. Coulter Department of Biomedical Engineering,
Institute for Electronics and Nanotechnology
Georgia Institute of Technology

CREATING THE NEXT®

CURRENT HEART MONITORS

- Clinical
 - Standard (12 lead)
 - Stress test
- Ambulatory
 - Holter monitor (3-5 lead)
 - Patch devices (single lead)
- Fitness trackers
 - Wrist
 - Chest strap



© MAYO FOUNDATION FOR MEDICAL EDUCATION AND RESEARCH - ALL RIGHTS RESERVED.
<https://www.mayoclinic.org/>



Photo <https://urgentmednetwork.com>



[CC BY-SA-NC](https://creativecommons.org/licenses/by-sa/4.0/)



[CC BY-SA](https://creativecommons.org/licenses/by-sa/4.0/)



LIMITATIONS

Form factor



- Holter monitor still Gold-standard
- Wires and rigid components



<https://www.mayoclinic.org/>


Gel electrodes

- Irritate skin
- Prone to motion artifacts





Contact impedance

- Tight straps
- Sensitive to sweat
- Skin strain







ADVANCES IN AMBULATORY ECG





- Conformable
- Adhesive
- Wireless
- Low-profile
- Light-weight

ADVANCES IN AMBULATORY ECG




Au (200 nm)
Cr (5 nm)
Polyimide (1.2 μm)
Elastomer (500 μm)


ADVANCES IN AMBULATORY ECG



ADVANCES IN AMBULATORY ECG




BioRadio (commercial)




BioStamp (commercial)



ADVANCES IN AMBULATORY ECG






ADVANCES IN AMBULATORY ECG



- High power consumption
- (heavy device)
- Large device footprint
- Not breathable


- Lower power consumption
- Weight reduction
- Smaller footprint and profile
- Breathable
- Strain-isolating

KEY ENABLING FEATURES




Ecoflex 30

Silbione



KEY ENABLING FEATURES

- Resist in-plane strain




Beam cross section $I = \frac{bh^3}{12}$

$$\delta_{max} = \frac{PL^3}{48EI}$$


$$h_{min} = \frac{PL}{2Eb^3\delta_{max}}$$

$$\delta_{max} = \frac{5qL^4}{384EI}$$


- Remain flexible out-of-plane








$$h_{max} = \sqrt[3]{\frac{5qL^4}{32Ebr[1 - \cos(L/2r)]}}$$








KEY ENABLING FEATURES




KEY ENABLING FEATURES



KEY ENABLING FEATURES



REAL-TIME ALL-DAY MONITORING



SUMMARY


Traditional



Today



Next Generation



ACKNOWLEDGMENT



Bio-Interfaced Translational Nanoengineering Group
PI: Professor Woon-Hong Yeo



Thank You!

- **Speaker contact info**
 - Email: ysk@me.gatech.edu
 - Website: www.ysk.bio
- **Yeo Research Group**
 - Website: yeolab.gatech.edu/