

Development of a Solid Phase Extraction Method for Cadmium and Nickel Using Stearic Acid Modified Activated Carbon

Nakano Minami¹, Mai Furukawa¹, Ikki Tateishi², Hideyuki Katsumata¹, Satoshi Kaneco¹

¹Department of Applied Chemistry, Graduate School of Engineering, Mie University, Japan

²Center for Global Environment Education & Research, Mie University, Japan

INTRODUCTION


Heavy metal

Heavy metals are widely used in various industrial applications.

Heavy metals have adverse effects through wastewater and emissions from factories.


Major uses of heavy metal

- Rechargeable batteries
- Automobiles
- Large machinery



Effects of heavy metal pollution

- Environmental pollution
- Health problems
- Ecosystem imbalance



Graphite Furnace Atomic Absorption Spectrophotometer (GFAAS)

- High sensitivity for measurement
- Easy operation
- Low cost of equipment

- Sensitive to physical and chemical interference
- Narrow quantitative range

Preconcentration of analytes prior to determination

Preconcentration

Purpose: To concentrate the analyte for measurement

Solid-phase extraction (SPE): A sample solution is passed through a solid phase, allowing separate the target material to be separated from impurities.

Purpose of this study

Preconcentration of trace cadmium and nickel using stearic acid modified activated carbon (SA-CAA)

METHOD

Adsorption preparation




Fig. 1. Preparation of SA-CAA.

Experimental

Solid phase extraction




Fig. 2. Solid phase extraction method.

RESULTS

Characterization

BET

Table 1. BET surface area, total pore volume, and average pore diameter of CAA and SA-CAA.

	V_m (cm^3 (STP) g^{-1})	a_s BET (m^2 g^{-1})	V_p (cm^3 g^{-1})	d_p (nm)
CAA	262.9	1140	0.96	3.3
SA-CAA	268.9	1170	0.97	3.3

V_m : Monolayer adsorption volume, a_s BET: BET Specific surface area, V_p : Full pore volume, d_p : Average maximum diameter.

SEM




Fig. 3. SEM images of (a) CAA and (b) SA-CAA.

RESULTS

Adsorbent type




Fig. 4. Examination of adsorbent type. Experimental conditions: stirring time 5 min, adsorbent amount 50 mg, initial pH 3, desorbent $0.1 \text{ mol L}^{-1} \text{ HNO}_3$ 3 mL, sample volume 200 mL

SA amount




Fig. 5. Examination of SA amount. Experimental conditions: stirring time 5 min, adsorbent amount 50 mg, initial pH 3, desorbent $0.1 \text{ mol L}^{-1} \text{ HNO}_3$ 3 mL, sample volume 200 mL

Stirring time




Fig. 6. Examination of stirring time. Experimental conditions: adsorbent amount 50 mg, initial pH 3, desorbent $0.1 \text{ mol L}^{-1} \text{ HNO}_3$ 3 mL, sample volume 200 mL

Adsorbent amount




Fig. 7. Examination of adsorbent amount. Experimental conditions: stirring time 20 min, initial pH 3, desorbent $0.1 \text{ mol L}^{-1} \text{ HNO}_3$ 3 mL, sample volume 200 mL

pH




Fig. 8. Examination of initial pH. Experimental conditions: stirring time 20 min, adsorbent amount 50 mg, desorbent $0.1 \text{ mol L}^{-1} \text{ HNO}_3$ 3 mL, sample volume 200 mL

Sample volume




Fig. 9. Examination of sample volume. Experimental conditions: stirring time 20 min, adsorbent amount 50 mg, initial pH 3, desorbent $0.1 \text{ mol L}^{-1} \text{ HNO}_3$ 4 mL

Table 2. Optimum conditions.

SA amount (CAA 1 g)	0.2 g
Stirring time	20 min
Adsorbent amount	50 mg
Initial pH	3
Desorbent	0.1 M HNO_3 4 mL
Sample volume	300 mL

Table 3. Analytical performance.

	Cd	Ni
Linear range (ppb)	0.01~0.1	0.1~1
Correlation coefficient (R^2)	0.99	0.99
Limit of detection ($3S^{(a)}/N^{(b)}$, ppt)	3	205
Limit of quantitation ($10S/N$, ppt)	21	684
RSD^(c) (%) (Cd:0.05 ppb, Ni:1 ppb, n=9)	12.6	14.3

(a) signal, (b) noise, (c) relative standard deviation.

CONCLUSION

- ✓ The use of stearic acid-modified activated carbon as the adsorbent showed selectivity for Cd and Ni, with improved recovery compared to unmodified activated carbon.
- ✓ The detection limits were 3 ppt for Cd and 205 ppt for Ni in this method.

FUTURE WORK / REFERENCES

- ✓ The effects of matrix elements (e.g., Mg, Ca, Na, K) on this method will be investigated, and its application to real samples will be studied.