



Universidad de Oviedo

GIVING NEW USES TO GLUCOSE METERS: DETECTION OF PROSTATE CANCER

IECB
2020

CLARA ABARDÍA-SERRANO

R. Miranda-Castro, N. de-los-Santos-Álvarez,

M. Jesús Lobo-Castañón

Emerging techniques

LIQUID BIOPSY

- Potential supplement of a traditional biopsy for early cancer diagnosis
- Monitoring the level of tumor biomarkers present in body fluids



Point-of-care (POC) Testing,

“Medical diagnostic test carried out outside the clinical laboratory, close to the place where the patient is being treated”

- Immediate, portable, easy to use and low-cost analyses



Pregnancy test

- Human chorionic gonadotropin detection
- Semiquantitative results



Glucose meter

- Glucose detection
- Quantitative results

Personal glucose meter (PGM)

- Most commonly known and widely developed point-of-care testing
- Device to control blood glucose concentration in people with diabetes
- Electrochemical detection: oxidoreductase enzymes

- Promising portable meter for targets beyond glucose
- The aim of this challenge: to generate glucose as a final product
 - Encapsulation of glucose solutions in nanocontainers
 - Enzyme-catalyzed glucose production



New targets

- Prostate cancer
 - Most common cancer in men
 - Damages prostate tissue
 - The probability of getting it increases with age

- 🎯 Prostate Cancer Antigen 3 or PCA3
 - 🎯 Urinary nucleic acid biomarker for prostate cancer diagnosis
 - 🎯 Synthetic DNA target of PCA3 RNA




Objective


- The main objective is to transform **personal glucose meter** into an electrochemical transducer of hybridization event, by introducing alkaline phosphatase onto magnetic beads as a **glucose-generating enzyme**, for the detection of the urinary biomarker for prostate cancer PCA3.

Proof of concept assay: $A_{20}-T_{20}$


1. Magnetic beads modification: Sensing layer construction



2. Hybridization genoassay onto magnetic beads






Dynamic range:
20 to 600 mg glucose dL^{-1}



Trinder method

- Enzymatic method for the determination of glucose
- Main enzyme: glucose oxidase, Gox
- Enzyme for colorimetric determination: horseradish peroxidase, HRP
- Spectrophotometric recorded, maximum absorbance at 505 nm



Optimization



Biotin - streptavidin

Target Poli A	50 nM
Hydrolysis buffer*	0.1 M Na ₂ CO ₃ de pH 9.9
Shaking*	1300 rpm
Temperature*	37°C
Glucose-1-phosphate*	3 g L ⁻¹
Hydrolysis time*	2.5 h
Strep-ALP	20 mg L ⁻¹

Glucose-6-phosphate	DEA	129 ± 11 mg dL ⁻¹
	Na ₂ CO ₃	130 ± 5 mg dL ⁻¹
Glucose-1-phosphate	DEA	163 ± 4 mg dL ⁻¹
	Na ₂ CO ₃	149 ± 3 mg dL ⁻¹

*DEA: buffer of 0.1 M of diethanolamine, pH 9.8, 1 mM MgCl₂


*J. Diabetes Sci. Technol. 2014, 8 (4), 855-858



Optimization

Biotin - streptavidin


Target Poli A	50 nM
Hydrolysis buffer	DEA
Shaking*	1300 rpm
Temperature*	37°C
Glucose-1-phosphate	30 g L ⁻¹
Hydrolysis time	2 h
Strep-ALP	20 mg L ⁻¹



Effect of the multivalence of streptavidin?

Stoichiometry of the enzymatic conjugate 2:1

(2 streptavidin : 1 molecule of alkaline phosphate)





Genoassay for PCA3



Genoassay

antifluorescein - Fab

1. Magnetic beads modification: Sensing layer




2. Sandwich-type genoassay: Homogeneous hybridization





Genoassay for PCA3


3. Sandwich-type genoassay: Heterogeneous hybridization



Results

antifluorescein - Fab

- ✓ Linear range: 10 – 100 pM
- ✓ Sensitivity (mg glucose dL⁻¹pM⁻¹):
 - ✓  PGM: 2.8 ± 0.1
 - ✓  Trinder method: 1.81 ± 0.05
- ✓ LOQ: 5 femtomoles in 0.5 mL of sample





Amplification strategy



Amplification

antifluorescein - Fab

1. Magnetic beads modification: Sensing layer



2. Sandwich-type genoassay: Homogeneous hybridization




Amplification strategy

Amplification

antifluorescein - Fab


3. Sandwich-type genoassay: Heterogeneous hybridization



Results

antifluorescein - Fab

- ✓ Linear range: 5 – 100 pM
- ✓ Sensitivity (mg glucose dL⁻¹pM⁻¹):
 - ✓ Two signaling probes: 3.3 ± 0.2
 - ✓ One signaling probe: 2.8 ± 0.1
- ✓ LOQ: 2.5 femtomoles in 0.5 mL of sample



- ✓ Selectivity assay against the PSA interferent

✓ 100 pM PSA interferent = $<20 \text{ mg glucose dL}^{-1}$ → Excellent selectivity

✓ 100 pM PCA3 = $343 \pm 35 \text{ mg glucose dL}^{-1}$

Conclusions

- A sandwich genoassay for quantification of PCA3 has been developed by introducing alkaline phosphatase as a tracer, through a label system on a tagged hybrid formed onto magnetic beads.
- Picomolar concentrations of PCA3 were amplified into millimolar concentrations of glucose by a monovalent label system, that can be detected with the personal glucose meter.
- The use of two signaling probes allowed the incorporation of two ALP molecules per analyte, resulting in an experimental LOQ of 5 pM.
- Excellent selectivity against a DNA homolog of PSA mRNA, also present in urine of prostate cancer patients, is achieved.

Acknowledgements

Electroanalysis Research Group, University of Oviedo

MSc. Raquel Sánchez Salcedo

MSc. Ramón Lorenzo Gómez

MSc. Ana Díaz Fernández

MSc. Clara Abardía Serrano

Dr. Rebeca Miranda Castro

Dr. Noemí de los Santos Álvarez

Dr. María Jesús Lobo Castañón





Universidad de Oviedo

GIVING NEW USES TO GLUCOSE METERS: DETECTION OF PROSTATE CANCER

IECB
2020

CLARA ABARDÍA-SERRANO

R. Miranda-Castro, N. de-los-Santos-Álvarez,

M. Jesús Lobo-Castañón