



Abstract Purification of Type II Collagen from Chicken Sternum *

Cheima Djehiche *, Nadia Benzidane and Lekhmici Arrar

Laboratory of Applied Biochemistry, Faculty of Nature and Life Sciences, Ferhat Abbas University, Setif 1., Setif 19000, Algeria; nadia19000@hotmail.fr (N.B.); lekharrar@univ-setif.dz (L.A.)

- * Correspondence: djehichecheima@gmail.com
- Presented at the 1st International Electronic Conference on Processes: Processes System Innovation, 17–31 May 2022; Available online: https://ecp2022.sciforum.net.

Abstract: The study presented in this work focused on the extraction and purification of type II collagen from chicken sternum using an acid-enzymatic extraction. Cartilage digestion with pepsin and acetic acid was followed by several successive precipitations using NaCl. Collagen was then solubilized and recovered by acetic acid. The estimation of the purity of obtained protein was carried out by polyacrylamide gel electrophoresis in presence of SDS and UV-visible scanning. Results of SDS-PAGE showed two pure proteins corresponding to α and β chains of collagen type II. Spectrophotometric scan gives a characteristic spectrum of collagen II with a single peak at 225 nm in absence of other proteins which absorb at 280 nm. The purified collagen from chicken sternum cartilage could have several laboratory applications especially in the induction of rheumatoid arthritis.

Keywords: collagen type II; chicken sternum; pepsin; precipitation; electrophoresis; absorption spectrum

Citation: Djehiche, C.; Benzidane, N.; Arrar, L. Purification of Type II Collagen from Chicken Sternum. *Proceedings* **2022**, *69*, x. https://doi.org/10.3390/xxxx

Academic Editor(s):

Published: date

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2022 by the authors. Submitted for possible open access publication under the terms and conditions of the Creative Commons Attribution (CC BY) license (https://creativecommons.org/licenses/by/4.0/).