

Universidade do Minho

# **Evaluation of essential oil-loaded** nanofibrous mats against the Escherichia virus MS2, a mimic of SARS-CoV-2, for prospective personal protective equipment uses





# CENTRO DE CIÊNCIA E **TECNOLOGIA TÊXTIL**

www.2c2t.uminho.pt

# S.F. Silva, J. Domingues, J.C. Antunes and Helena P. Felgueiras\*

Centro de Ciência e Tecnologia Têxtil (2C2T), Universidade do Minho, Portugal \*helena.felgueiras@2c2t.uminho.pt

### Introduction

In December 2019, a novel strain of coronavirus, SARS-CoV-2, was identified. Infected patients revealed symptoms of fever, cough (dry), sore throat, and fatigue, which began manifesting after 5 days of incubation. Hoping to prevent transmission, many countries adopted a mandatory mask use in closed public spaces. However, most mask options display a passive action against COVID-19. To overcome such restrictions, this work proposes the incorporation of anti-viral essential oils (EOs) loaded onto a nanofibrous layer that can be adapted to both community and commercial masks.


### **EOs Drawbacks**

- cytotoxic at increased concentrations, which prevents systemic delivery;
- present low resistance to degradation by external factors (e.g., temperature, light, moisture);
- highly volatile in their free, unloaded form.


# **Results and Discussion**

# Nanofiber Loading

### Strategy 1 Physical adsorption of the EOs at the surface of the electrospun PCL mats (label PCLaEOs)



### Nanofiber Morphology




## **Materials and Methods**

Syringe or Capillary Tube

### **Electrospinning**


Spinning technique that allows the production of continuous, homogeneous nanofibers.



Collecting

Plate

# **Polymeric solution preparation** Polycaprolactone (PCL)



### Why PCL?

- Food and Drug Administration (FDA)-approved;
- Drug loading efficiency (i.e., EOs biomolecules);
- Excellent mechanical strength, non-toxic, hydrophobic, biodegradable, etc.

PCL at 14 wt.% in chloroform/dimethyl formamide (CHF/DMF at 9/1 v/v)


### **Electrospinning processing conditions**

Potential: 23 kV Extruding Speed: 0.7 mL/h Distance to Collector: 26 cm Needle (inner diameter): 18 gauge

### **Essential Oils (EOs) Selection**

20 EOs\* with antimicrobial potential were examined for their minimum inhibitory concentrations (MICs) against the MS2 Escherichia host and for

### **Confirmation of EOs Loading**



## **Antimicrobial Effect**

their virucidal concentration (VC) against the MS2 virus, mimic of SARS-Cov-2, at initial concentration of 1x107 CFUs or PFUs/mL, respectively.

EOs	MICs against MS2 host (mg/mL)	VC against MS2 virus (mg/mL)
Lemongrass (LGO)	178.0	356.0
Niaouli (NO)	45.7	365.2
Eucalyptus (ELO)	>	<b>586.0</b>
Orchid	85.6	428.0
Tea Tree Oil	22.4	447.5
Clove	105.6	528.0

\*DOI: 10.3390/biom10081129


### Acknowledgments

This work is financed by FEDER funds through COMPETE and by national funds through FCT via the projects POCI-01-0145-FEDER-028074 and UID/CTM/00264/2021.

**Contact Killing Examinations** Microorganism: Bacteriophage MS2 (mimic of SARS-CoV-2) Initial Concentration (t = 0 h):  $3.8 \times 10^7$  PFUs/mL Incubation Period: 4 h Temperature: room

MS2 reduction was observed after 4 h of interaction on all EOloaded mats. PCL blended with LGO was the most effective from the group in fighting the virus.

PCL on its own was also seen to retain virus within its structure, attesting to its functionality as a potential retaining layer for masks.



**Conclusions:** Data demonstrated the potential of these EOs-loaded PCL/CA nanofiber mats to work as COVID-19 active barriers for individual protection masks.

