Genetically encoded fluorescent probes for imaging of intracellular localization and activity of SARS-CoV-2 proteins

This work was supported by the Russian Foundation for Basic Research grant 20-04-60370

Elena L. Sokolinskaya, L.V. Putlyaeva, K.A. Lukyanov

Center of Life Sciences, Skolkovo Institute of Science and Technology, Moscow, Russia

Skoltech

Results

elena.sokolinskaya@gmail.com



1.3-1.6-fold increase in the mScarlet/iRFP670 signal ratio several hours after induction of

PLpro expression



Since December 2019, the problem caused by severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) has grown to a global threat. The search for new treatment strategies is strongly associated with both fundamental research on mechanisms of virus life cycle and development of new screening platforms for anti-viral drug candidates.

In this project we labelled SARS-CoV-2 membrane proteins M, E, S to track protein intracellular localization and study protein-protein interactions in mammalian cells using fluorescent microscopy.

In addition, we developed several types of FRET-based and translocation sensors to track SARS-CoV-2 PLpro protease and measure its activity in living cells.


The proposed strategy for studying viral proteins combines two important advantages. Firstly, research is conducted on human living cells, which is closely approximated to native conditions in contrast to in vitro experiments. Secondly, the experimental system lacks interaction with a functional virus which makes it completely safe for the researcher.

Conclusions


1. E protein accumulates in cell lysosomes and endosomes. S protein is distributed in the form of a network (obviously, endoplasmic reticulum).


2. M protein localizes in two compartments – Golgi and lysosomes. Untagged M shows a tendency for oligomerization.

3. Two types of sensors (soluble FRET-sensor and membraneassociated translocation sensors) for real-time detection of PLpro activity in live cells are developed. Testing in mammalian cells demonstrated a robust readout.



Translocation sensor for SARS-CoV-2 PLpro activity





I) changes its intracellular localization: the ER network disappears and the red signal becomes distributed in the cytoplasm. HeLa cells.