

Universidade do Minho



Cinnamon leaf oil release from chitosan microcapsules embedded within a sodium alginategelatin hydrogel-like film for the inhibition of multidrug resistant bacteria

Catarina S. Miranda*, Joana C. Antunes, Natália C. Homem, Helena P. Felgueiras.

Centro de Ciência e Tecnologia Têxtil (2C2T), Universidade do Minho, Portugal *catarina.miranda@2c2t.uminho.pt

CENTRO DE CIÊNCIA E **TECNOLOGIA TÊXTIL**

www.2c2t.uminho.pt

Introduction

Pseudomonas aeruginosa is considered a public threat, with a great number of antibiotic resistant microorganisms. Essential oils (EOs) have demonstrated significant effects against several pathogens, however, they present a volatile nature, therefore, biodegradable polymeric delivery platforms are being engineered. Here, hydrogel-like films were produced from a combination of sodium alginate (SA) and gelatin (GN) to serve as delivery platforms for the controlled release of cinnamon leaf oil (CLO) entrapped within chitosan (CS) microcapsules. CS microcapsules were prepared via ionotropic gelation with tripolyphosphate (TPP), with CLO at minimum inhibitory concentration (MIC) in the core. Microcapsules were then embedded within a biodegradable SA/GN polymeric matrix processed via a solvent casting/phase inversion methodology with SA/GN used at 70/30 polymer ratio and 2 wt% SA concentration in distilled water. The coagulation bath was composed of a 2 wt% calcium chloride (CaCl₂) aqueous solution. Qualitative and quantitative antimicrobial examinations validated the modified film potential to fight infections caused by *P. aeruginosa* bacteria.

Goal

Engineer a hydrogel-like film delivery platform that guarantees the controlled release of cinnamon leaf oil for the eradication of *P. aeruginosa*-derived infections.

CLO Encapsulation: Chitosan microcapsules



Films' morphology



Very hydrated films; Flexible; Homogeneous; Regular distribution of microcapsules

Chemical characterization - FTIR



Sample	CS (mg)	dH₂O (mL)	Acetic acid (µL)	CLO (µL)	TPP (mL)	pH adjust	NaOH (mL)
CS1	2.5	9.9	100	834	9	No	-
CS4	2.5	8.8	100	834	9	Yes	1.1





CS1/CS4 solutions with **Microcapsules**

Ultraturrax (IKA T25) 10000 rpm, 10", 5 x



Purification by

dialysis



Final CS1/CS4 solutions with Microcapsules --- SA/GN's characteristic peaks

CLO release profile



inferior throughout the test on the **CS4-loaded films**

CLO release

profile was

Antimicrobial action



Enhanced performance of the SA/GN/CS1

Films Production



The incorporation of the CLO-containing CS microcapsules within the SA/GN fibers was confirmed

The continuous release of the entrapped oil over a period of 24h was attained, with a matched time kill kinetics against *P. aeruginosa*

CS1 loaded films were determined more effective than CS4 loaded films

Acknowledgments

This work is financed by FEDER funds through COMPETE and by national funds through FCT via the projects POCI-01-0145-FEDER-028074 and UID/CTM/00264/2021. CSM also acknowledges FCT for PhD grant (2020.08547.BD).

